



TEST REPORT

Reference No. : WTF20F12095227C
Applicant : Mid Ocean Brands B.V.
Address : 7/F., Kings Tower, 111 King Lam Street, Cheung Sha Wan, Kowloon, Hong Kong
Manufacturer..... : 115253
Sample Name..... : Double wall vacuum bottle
Model No...... : MO9431
Test Requested..... : Client's requirements - Dishwasher safe test.
Test Method : Please refer to next page (s)
Test Conclusion : Please refer to next page (s)
Date of Receipt sample..... : 2020-12-10
Date of Test..... : 2020-12-10 to 2020-12-15
Date of Issue..... : 2020-12-16
Test Result..... : Please refer to next page (s)
Note : As specified by client, only test the designated sample.

Remarks:

The results shown in this test report refer only to the sample(s) tested; this test report cannot be reproduced, except in full, without prior written permission of the company. The report would be invalid without specific stamp of test institute and the signatures of compiler and approver. If the report is not stamped with the accreditation recognized seal, it will only be used for scientific research, education, and internal quality control activities, and is not used for the purpose of issuing supporting data to the society.

Prepared By:

Waltek Testing Group (Foshan) Co., Ltd.

Address: No.13-19, 2/F., 2nd Building, Sunlink International Machinery City,
Chencun, Shunde District, Foshan, Guangdong, China

Tel:+86-757-23811398 Fax:+86-757-23811381 E-mail:info@waltek.com.cn

Compiled by:

Rena Chen

Rena.Chen / Project Engineer

Approved by:

Swing Liang

Swing.Liang / Technical Manager

**Test Result:****Dishwasher safe test (In House Method)**

No. of Specimen	:	1 piece for testing. Remains as control sample
No. of Cycles	:	(10) cycles
Dishwasher	:	Siemens SD6P1S
Detergent Brand	:	Cascade
Washing condition	:	
1. Cycle Time	:	~ 90 Minutes.
2. Cycle Water Temperature	:	Main Wash ~ 140°F (60°C) Rinse ~ 161°F (72°C)
3. Location inside dishwasher	:	Bottom rack
Testing Procedure	:	1) Visually inspect sample and record any aesthetic defects. 2) Load into the dishwasher and run the wash cycle as per the above condition. 3) Examine from a distance of ~ 30cm and record any difference as compared with the control sample. 4) Repeat 2) ~ 3) until the required number of cycle is completed.
Test Requirement	:	No discoloration, deformation, staining of water (water mark), or substantial deterioration of item's appearance and/or function.
Conclusion(s)	:	PASS 1. No noticeable change in appearance / size was found on the sample(s) 2. The tested sample(s) remained functional

Test Specimen Description:

No.1: Double wall vacuum bottle (Overall)

WALTEK





Sample photo:



==== End of Report =====

WALTEK

